nourish seating | product specification

Healthcare Patient Chairs

Models: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20, GC3685,

GC3698

Healthcare Ottoman Model: GC3686

Dimensions Standard Chairs

(Models: GC3681, GC3682, GC3683, GC3684, GC3698)

•		•
Overall Depth		29 in / 737 mm
Overall Width		28.5 in / 724 mm
Overall Height		45.5 in / 1156 mm
Seat Depth		18.5 in / 470 mm
Seat Width		20 in / 508 mm
Seat Height	GC3681, GC3682, GC3683 GC3698	17-20 in / 432-508 mm 17 in / 432 mm
Back Height - From Seat		29 in / 737 mm
Armrest Clearance		21 in / 533 mm
Arm Height - From Seat		9 in / 229 mm
Product Weight		42 lb / 19.1 kg
Shipping Cube		33.5 ft ³ / 0.95 m ³

High Seat Chairs

(Models: GC3681-20, GC3682-20, GC3683-20, GC3684-20)

Overall Depth	29 in / 737 mm
Overall Width	28.5 in / 724 mm
Overall Height	45.5 in / 1156 mm
Seat Depth	18.5 in / 470 mm
Seat Width	20 in / 508 mm
Seat Height	20-23 in / 508-584 mm
Back Height - From Seat	29 in / 737 mm
Armrest Clearance	21 in / 533 mm
Arm Height - From Seat	9 in / 229 mm
Product Weight	46 lb / 20.9 kg
Shipping Cube	33.5 ft ³ / 0.95 m ³

Bariatric Chair (Model: GC3685)

Overall Depth	29 in / 737 mm
Overall Width	40 in / 1016 mm
Overall Height	47 in / 1194 mm
Seat Depth	18.5 in / 470 mm
Seat Width	35 in / 889 mm
Seat Height	18.5 in / 470 mm
Back Height - From Seat	29 in / 737 mm
Armrest Clearance	21 in / 533 mm
Arm Height - From Seat	9 in / 229 mm
Product Weight	60 lb / 27.3 kg
Shipping Cube	45 ft³ / 1.27 m³

Ottoman (Model: GC3686)

Overall Depth	17.5 in / 445 mm
Overall Width	17 in / 432 mm
Overall Height	12.5 in / 318 mm
Product Weight	15 lb / 6.8 kg
Shipping Cube	3.6 ft ³ / 0.10 m ³

Contact sales representative for more information.

GC3686 GC3681-20 GC3685





nourish seating | product specification

Features

- · Transfer arms fold away to allow easy access for transfer of less mobile users
- Easy grip, forward extending Self-Skinned Urethane (SSU) armcaps

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20, GC3685

• Chairs have casters on rear legs to assist chair movement

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20

• Height adjustable legs provide 3 in (76 mm) seat height range or seat angle adjustment

MODELS: GC3682, GC3682-20, GC3684, GC3684-20

• Schukra mechanism provides lumbar depth adjustment

MODELS: GC3683, GC3683-20, GC3684, GC3684-20

• Schukra mechanism integrated in headrest provides depth adjustment

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20

• Two rear tip and move wheels allow easy relocation of unoccupied chairs

MODEL: GC3698

• Four wheel chair with fixed length legs

MODEL: GC3686

• Matching ottoman

Options

- Frame is standard in Tungsten. Also available in six other colors
- Armcaps are standard in Black. Also available in five other colors, matching frame

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20

• Easy grip bar can be added to assist chair movement. Available in four colors

COM Yardage

Yardage is based on unmatched materials with widths of 54 inches for seating

MODELS: GC3681,

GC3681-20, GC3682,

GC3682-20, GC3683,

GC3683-20, GC3684,

GC3684-20, GC3698

• Back: 1.30 yards

• Seat: 0.70 yards

• Total: 2.00 yards

MODEL: GC3685

• Back: 1.00 yards

• Seat: 2.25 yards

• Total: 3.25 yards

MODEL: GC3698

• Back: 1.30 yards

• Seat: 0.70 yards

• Total: 2.00 yards

Infection Control

- · No Velcro on seat or back
- Seam sealing option
- Moisture barrier option
- Integrated Schukra mechanisms (in applicable models) allow back and headrest adjustment with an easy to clean single piece back
- Full 360° cleanout, 0.75 in (19 mm) wide
- · Back is removable and field replaceable
- For proper cleaning and disinfecting practices, the textile selection needs to be identified. It is vital that following any chemical application a clear water wipe down occur to avoid chemical build up that will lead to deterioration of the textile. Contact Global should information on using a specific disinfectant be required

Construction (Overall Frame)

- 3 in (76 mm) diameter steel tube legs, 11 gauge
- All metal-to-metal connections

Construction (Seat)

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20, GC3698

• 0.5 in (13 mm) thick plywood seat pan, 7-ply

MODEL: GC3685

• 0.625 in (16 mm) thick plywood seat pan, 10-ply

Construction (Back)

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20, GC3698

- 0.5 in (13 mm) thick plywood back, 6-ply
- · Synthetic suspension webbing

MODEL: GC3685

• 0.375 in (10 mm) thick plywood back pan, 7-ply

Construction (Arms)

- · Tungsten color transfer arms
- SSU armcaps

Construction (Glides)

- Non-marking black glides
- Injection molded polypropylene glides

Construction (Casters)

• 3 in (76 mm) diameter caster

Upholstery

· Seat and back can be upholstered in fabric, vinyl or leather

Foam (Seat)

- Dual layer, Ultracell (high density) seat foam
 - 1.5 in (38 mm) thick layer, Density: 2.5 lb. ILD: 50
 - 1 in (25 mm) thick layer, Density: 2.6 lb. ILD: 35





nourish seating | product specification

Foam (Back)

- 1.2 in (30 mm) thick back foam
- Density: 1.9 lb. ILD: 39

Weight Capacity

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20

- Designed for an active weight capacity of 350 lb MODEL GC3698
- Designed for an active weight capacity of 500 lb MODEL GC3685
- Bariatric chairs have active weight capacity of 750 lb

Warranty

MODELS: GC3681, GC3681-20, GC3682, GC3682-20, GC3683, GC3683-20, GC3684, GC3684-20, GC3698

Lifetime warranty on all components except the following:

- 5 year warranty on foam and textiles
- 12 year warranty on control mechanisms
- · Warranty applies only to single shift, standard co mmercial usage (8 hours/day, 5 days/week)

MODEL: GC3685

12 year warranty on all components for 24 hour use except the following:

• 5 year warranty on foam and textiles

Fire Rating

All textiles and foam is tested and certified, meeting California Technical Bulletin 117 (CAL 117-2013) standards

Environmental

View Environmental Data Sheets online.

All foam is produced without using any ozone depleting chloro-flouro chemicals.

Awards/Certifications

- Meets or exceeds all applicable BIFMA standards for safety and durability
- Manufactured in ISO 9001 and 14001 certified facilities

Please visit globalfurnituregroup.com for additional product information including environmental certifications.

